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9. Sources/Quellen

- Figiel L, Kasparzek JD, Peruga J, Lipiec P, Drozd P, Krzeminska-Pakula M, Smigielski J:** Heart-fatty acid binding protein – a reliable marker of myocardial necrosis in a heterogeneous group of patients with acute coronary syndrome without persistent ST elevation. *Kardiol Pol* (2008), 66:253-259.
- Adams, et al.** Diagnosis of Perioperative myocardial infarction with measurements for cardiac troponin I. *N.Eng.J.Med.* 1994, 330:670.
- Mensink, M et al;** Lifestyle changes and lipid metabolism gene expression and protein content in skeletal muscle of subjects with impaired glucose tolerance. *Diabetologia* 2003, 46: 1082.
- Alpert JS, et al.** Myocardial Infarction Redefined, Joint European Society of Cardiology American College of Cardiology: *J. Am. Coll. Cardio.*, 36(3)959, 2000
- Hossein-Nia M, et al.** Cardiac troponin I release in heart transplantation. *Ann. Thorac. Surg.* 1996, 61: 227.
- Cubranic Z. et al.,** Diagnostic accuracy of heart fatty acid binding protein (H-FABP) and glycogen phosphorylase isoenzyme BB (GPBB) in diagnosis of acute myocardial infarction in patients with acute coronary syndrome, *Biochemia Medica* 2012; 22(2): 225–36.

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
Product for single-use only
Einmalprodukt

Note: Read the operating instructions of the product before each use of the test. The kit is for medical professionals only.

4. Evaluation

A heart attack is present if the figure below shows the following:

Positive Event:

In the matrix displayed above you can see the potential results of the BelCARDIO test and their correlation.

- ❑ A positive hFABP value indicates a heart attack more than 20 minutes and 24 hours after the onset of the clinical symptoms.
- ❑ If hFABP and Troponin I are both positive, an infarct has a potential occurrence within the last 4 to 24 hours.
- ❑ If only Troponin I is positive, it is highly probable that a heart attack has occurred within the last 24 hours to 8 days.

Negative Event:

A pink coloured line appears only in the control line region (C). The test has performed the necessary functions.

In the two test line areas (cTnI) and (hFABP) no recognizable pink coloured line is present. The test detected neither Troponin I nor hFABP.

5. Performance Characteristics

	hFABP	cTroponin I
Rel. Sensitivity	89,9%	99,4%
Rel. Specificity	91,0%	99,0%
Accuracy	90,7%	99,1%

1. Introductory Remarks

Weaknesses in infarct diagnostics

Troponin I can be detected in the blood only approx. 3 hours after an ischemic event. Other biomarkers (e.g. myoglobin) are not sufficiently reliable because they are not specific to the body. Even a 12-channel ECG remains without any detection in about 25% of infarct patients. In some cases, an acute infarct could not be clearly diagnosed.

Concept of BelCARDIO


In order to improve the sensitivity in acute infarcts and weak NSTEMI, BelCARDIO determines the early phase protein hFABP as well as troponin I. Due to its biochemical properties and physiological function, the protein hFABP is particularly well suited to the detection of myocardial infarction within the first 12 hours after the onset of symptoms. Only one sample introduction allows the largest possible diagnostic time frame. Few minutes for your diagnosis give a double security, thus the chances of survival of your patient rise clearly.

Clinical Performance

The hFABP, the heart specific fatty acid binding protein, is involved in the active fatty acid metabolism in which it transports fatty acids from the cell membrane of the mitochondria to the metabolic conversion. The clinical detection limit is 8.0 ng/mL.

The release occurs 20 minutes after the infarct. The increased concentration remains for about 24 hours and then falls back to the natural level.

Cardiac troponin I (cTnI) is a protein that occurs as a complex with troponin T and troponin C in the heart muscle. The components are released in the bloodstream and can be detected there several hours after the onset of a heart attack. Troponin I levels remain elevated for 8 days after release. Thus a larger time window results for the proof of a heart damage. The clinical detection limit is 0.5 ng/mL.

6. Test Accessories

Each BelCARDIO Test contains a cassette of pre-coated test strips with colour conjugates and reactive reagents. The sample field (S) is designed for the introduction of whole blood, capillary blood, plasma or serum sample (50µL).


Figure 1: New test cassette with test strips and markings.

The **detection strips** (hFABP) and (cTnI) on the lateral flow test are each coated with an antibody. The **control strip** (C) is coated with a specific polypeptide antibody (see Figure 1).

Accessories include a transfer pipette for applying the sample, an ampoule containing the buffer solution, a sterile lancing device for the fingertip and a user manual to complete the product. The color strips on the Lateral Flow Stripe are evaluated and displayed by a suitable scanner/reader at the correct time.

The lancing device must be unlocked and the fuse removed before use.


Figure 2: Prepare lancing device

Figure 3: Unlocking

Figure 4: Finger pricking

Figure 5: Transfer Blood

Figure 6: Pipette

Figure 7: Test Unit

2. Sample Preparation

BelCARDIO can be performed with whole blood from the fingertip, directly from a venipuncture site, plasma or serum.

Whole blood: If whole blood is not given directly from the fingertip or venipuncture site, it must be anticoagulated with citrate, EDTA or heparin. Whole blood samples from the vein should be stored at 2°-8°C when the test is performed within 2 days of sample collection. Whole blood samples must not be frozen.

Serum: Serum samples may be stored at 2°-8°C for up to three days. For longer storage, serum samples should be stored at -20°C or below. Bring samples to room temperature before testing.

Plasma: Separate the plasma as soon as possible to avoid hemolysis. Use only clear, non-haemolysed sample material. When using plasma samples, citrate, EDTA and heparin can be used as anticoagulants.

Moisture and incorrect temperature may adversely affect the test results.

Icteric, lipemic, hemolyzed, heat-treated and contaminated specimens may give false results. Repeat the test with a serum or plasma sample from the same patient using a new test device.

Do not store specimens at room temperature for extended periods of time. Plasma samples can be stored at 2°-8°C for up to three days.

For longer storage, plasma specimens should be stored at -20°C or below. Bring samples to room temperature before testing.

Frozen samples should be completely thawed and well mixed before testing. Samples should not be repeatedly frozen and thawed.

When samples are shipped, they should be packaged in accordance with all applicable regulations for the transport of etiologic agents.

There is a small possibility that some whole blood samples with a very high viscosity or whole blood samples that have been stored for more than 2 days may not be collected in the sample field of the test.

7. Test Properties

Storage and Stability: This test may be stored at 2 to 30 °C until the imprinted expiration date. Do not freeze the test. The test must remain in the sealed blister pack until use. Please protect the test components from contamination and moisture.

The test reagents and samples: Should be brought to room temperature (15-30°C) before use. After opening the aluminum foil bag, the tests should be performed as soon as possible.

8. Restrictions on Use

Application: This product is intended for in vitro diagnostic use under professional supervision. This test is a disposable product. Please handle it with care to avoid hazards and infections.

Validity and packaging: Keep product packaging intact before use. Please do not use products that exceed the validity period or are damaged. Moisture and incorrect temperature may adversely affect test results.

Warning: The test is intended only as an aid to diagnosis. The determination of the biomarkers cannot replace any of the applied clinical examination methods. The result is for clinical reference only. Treatment of patients should be based on a combination of symptoms, clinical signs, medical history, other laboratory tests and therapeutic responses.

Please follow the instructions carefully. Do not end the test halfway.

Precautions: Do not eat, drink, or smoke while handling specimens and performing the test. All specimens are considered potentially hazardous and must be treated as infectious specimens. Follow proven precautions against microbiological risks during the test procedure. Wear protective clothing such as disposable gloves and eye protection when testing samples.

Substances used: The test also includes products of animal origin. Certificates on the origin and/or health status of the animals cannot completely guarantee the absence of transmissible pathogenic components. It is therefore recommended to consider these products as potentially infectious due to the blood samples used and to take the appropriate precautions.

3. Test Execution

Each BelCARDIO Kit contains a cassette with test strips and reactive reagents. The sample field (S) is designed for the introduction of whole blood, capillary blood, plasma or serum sample (50µL).

1. Bring the patient's test, buffer solution and specimen material to room temperature (15°-30°C).
2. Prepare the stopwatch.
3. Do not open the pouch until everything is ready to perform the test.
4. Remove the test strip seal and mark the sample ID on the test strip.
5. Place the test strip on a level table.
6. Using the **transfer pipette**, take **2 drops** of whole blood, or serum or plasma sample, or 2 drops of blood from the fingertip (~50µL) and place in the application field (S). Alternatively, two hanging drops of blood from the fingertip can be placed directly into the application field.
7. **Add one drop of buffer solution.** Avoid formation of air bubbles in the application field.
8. **Start the stopwatch.**
9. Read the result exactly 10 minutes after the sample has been added.

or

10. When using an appropriate reader: Insert the test cassette into the reader, press the start button. The device stops automatically and displays the measurement results.
11. Under no circumstances should results be read after 15 minutes.
12. Note the results for the patient file.

Important:

After approx. 30 seconds you can observe the chromatography on the test strip in the cassette.

A red-pink color front forms which moves across the membrane.

In rare cases this process may be disturbed. Then add another drop of buffer solution.

9. Sources / Quellen

Figiel L, Kasparzek JD, Peruga J, Lipiec P, Drozd P, Krzeminska-Pakula M, Smigielski J: Heart-fatty acid binding protein – a reliable marker of myocardial necrosis in a heterogeneous group of patients with acute coronary syndrome without persistent ST elevation. *Kardiol Pol* (2008), 66:253-259.

Adams, et al. Diagnosis of Perioperative myocardial infarction with measurements for cardiac troponin I. *N.Eng.J.Med.* 1994, 330:670.

Mensink, M et al; Lifestyle changes and lipid metabolism gene expression and protein content in skeletal muscle of subjects with impaired glucose tolerance. *Diabetologia* 2003, 46: 1082.

Alpert JS, et al. Myocardial Infarction Redefined, Joint European Society of Cardiology American College of Cardiology: *J. Am. Coll. Cardio.*, 36(3)959, 2000

Hossein-Nia M, et al. Cardiac troponin I release in heart transplantation. *Ann. Thorac. Surg.* 1996, 61: 227.

Cubranic Z. et al., Diagnostic accuracy of heart fatty acid binding protein (H-FABP) and glycogen phosphorylase isoenzyme BB (GPBB) in diagnosis of acute myocardial infarction in patients with acute coronary syndrome, *Biochemia Medica* 2012; 22(2): 225–36.

10. Symbols / Symbole

REF Catalogue number Bestellnummer	Package size Anzahl der Tests pro Set 1-5
LOT Batch code Chargencode	Expiry date Verfallsdatum
IVD In Vitro Diagnostic Medical Device Für die In-vitro-Diagnostik	Storage conditions Lagerungstemperatur 2°C - 30°C
Manufacturer Hersteller	STERILE Sterile Product Steriles Produkt OX20 1TU, UK,c 0120 OX20 1TU UK,
CE CE marked as per Annex III of Directive 98/79/EC of the European Parliaments and Council of Oct. 27. 1998. CE-Kennzeichnung gemäß Anhang III der Richtlinie 98/79/EG des Europäischen Parliaments und des Rates vom 27. Oktober. 1998.	
Consult Instructions for Use Siehe Gebrauchsanweisung	Product for single-use only Einmalprodukt