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### Instructions for use

### for health care professionals

#### In-vitro-Diagnostik

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Product for single-use only  
Einmalprodukt

Manufacturer  
Hersteller

CE marked as per Annex III of Directive 98/79/EC of the European Parliaments and Council of Oct. 27. 1998.

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## 15. Sources/Quellen

1. Johansen, J. S., P. Christoffersen, S. Møller, P. A. Price, J. H. Henriksen, C. Garbarsch and F. Bendtsen (2000). "Serum YKL-40 is increased in patients with hepatic fibrosis." J Hepatol 32(6): 911-920.
2. Linlin Yan, Yongqiong Deng, Jiyuan Zhou, Hong Zhao, Guiqiang Wang (2017). "Serum YKL-40 as a biomarker for predicting significant fibrosis and advanced fibrosis in chronic hepatitis B patients with normal or mildly elevated ALT." Hepatol Int (2017) 11 (Suppl 1): S1-S1093.
3. Kamal, S. M., B. Turner, Q. He, J. Rasenack, L. Bianchi, A. Al Tawil, A. Nooman, M. Massoud, M. J. Koziel and N. H. Afdhal (2006). "Progression of fibrosis in hepatitis C with and without schistosomiasis: correlation with serum markers of fibrosis." Hepatology 43(4): 771-779.
4. Wang J, Gao F, Mo F, Hong X, Wang H, Zheng S, Lin B (2009). "Identification of for Liver cancer through integrative secretome and transcriptome analysis." Proteomics-Clinical Applications, 3: 541-551, 2009.
5. Hood, L., J. R. Heath, M. E. Phelps and B. Lin (2004). "Systems biology and new technologies enable predictive and preventative medicine." Science 306(5696): 640-643.
6. Huang, H., T. Wu, J. Mao, Y. Fang, J. Zhang, L. Wu, S. Zheng, B. Lin and H. Pan (2015). "CHI3L1 Is a Liver-Enriched, Noninvasive Biomarker That Can Be Used to Stage and Diagnose Substantial Hepatic Fibrosis." OMICS, a Journal of Integrative Biology 19(6): 339-345.
7. Kumagai, E., Y. Mano, S. Yoshio, H. Shoji, M. Sugiyama, M. Korenaga, T. Ishida, T. Arai, N. Itokawa, M. Atsukawa, H. Hyogo, K. Chayama, T. Ohashi, K. Ito, M. Yoneda, T. Kawaguchi, T. Torimura, Y. Nozaki, S. Watanabe, M. Mizokami and T. Kanto (2016). "Serum YKL-40 as a marker of liver fibrosis in patients with non-alcoholic fatty liver disease." Sci Rep 6: 35282.
8. Lin Wang, Yameng Sun, Hong You, Jidong Jia (2017). "The decrease of serum CHI3L1 levels after antiviral therapy can monitor reversal of liver fibrosis." Hepatol Int (2017) 11 (Suppl 1): S1-S1093.
9. Hakala, B. E., C. White and A. D. Recklies (1993). "Human cartilage gp-39, a major secretory product of articular chondrocytes and synovial cells, is a mammalian member of a chitinase protein family." J Biol Chem 268(34): 25803-25810.

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### For Healthcare Professionals

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
### In-vitro-Diagnostics

Rapid test to evaluate the disease-risk of the liver.

Qualitative detection of the biomarker YKL-40 (CHI3L1).  
Quantitative detection with additional reader.

*This rapid diagnostic test is intended for use in the early detection of potential liver diseases. In the daily use this test enables the medical doctor a quicker diagnosis and a precise determination of the various phases of liver function.*

## 6. Test Procedures

**NOTE:** Please read the instruction manuals of the product and the scanner/reader before using the test. The kit is intended for use by healthcare professionals only.

**Test reagents and samples** should be brought to ambient temperature (15-30°C) before use. After opening the aluminum foil bag, the assay should be performed as soon as possible.

**Application:** The test is only intended as an auxiliary diagnosis for liver fibrosis and cirrhosis. The determination of YKL-40 (CHI3L1) cannot replace any of the clinical examination methods used.

## 7. Quantitative Measurement

This product is only intended for the **quantitative measurement** of the analytical substance. Further laboratory tests are also recommended for confirmatory diagnosis of the disorder.

## 8. Scanner for quantitative Evaluation

- SkanFlexi Lateral Flow Multi-Reader of Skannex AS, Oslo, Norwegian ([www.skannex.com](http://www.skannex.com)).
- BelREADER of PharmACT AG, Berlin, Deutschland.

## 9. Test Characteristics

**Stability:** The test met all requirements for linearity, repeatability, negative reference agreement, specificity and sensitivity of the test kit 10 days after annealing at 37° C. The test was performed at a temperature of 10° C after annealing.

**Specificity:** QC samples (QM-Control) of AFP at 400ng/mL and C reactive protein (CRP) at 100µg/mL all reported values of <79ng/mL for multiple measurements.

**Accuracy:** The QC sample (quality control) of CHI3L1 at 80ng/mL and 200ng/mL both have a detection range within 20% of the QC values.

**Linearity of dose-response curve:** Linearity range is between 25 and 800ng/mL. The correlation coefficient (R) of the fitted line of ratios and protocol concentrations is > 0,95.

**Reproducibility:** The coefficient of variation is less than 15% when the CHI3L1 QC is measured 10 times at 200ng/mL.

**Negative reference agreement:** 15 negative reference samples all reported values of <79ng/mL. Detection limit: 25ng/mL.

## Introduction

### Weaknesses in liver function testing

The liver takes an important part in the metabolism of the body, it is the main "detoxifying organ". Every day toxic substances are produced in the metabolism, which the liver converts into harmless substances and passes on to the body for excretion. In the diagnostics liver diseases are not immediately recognizable, only a few humans show clear symptoms, cause the liver functions usually by way of relatively well.

Liver disorders are extremely expensive in terms of visits to doctors and hospitals, premature loss of productivity and human suffering. In the diagnostics the liver is judged typically on the basis the determination of the blood values and, where unclear findings are present, by special ultrasonic investigations (e.g. Fibroscan).

### Risk factors for liver disorders

Excessive, high-fat foods combined with high alcohol consumption and lack of movement can lead to fatty liver. Diabetes mellitus can also promote this disorder. Approximately 40% of all liver diseases are induced by alcohol and/or overweight, including lack of exercise.

Hepatitis A and overconsumption of alcohol are the biggest risk factors for liver disease. Already a daily intake of 10g alcohol for women and 20g alcohol for men is a critical factor.

### Typical Liver Diseases

#### Non-alcoholic fatty liver disease (NAFLD)

Humans with a type 2 diabetes have a greater risk to suffering from severe liver diseases. The risk is 3x higher in men and 5x higher in women. NAFLD is often associated with adiposity, which is also a risk factor for type 2 diabetes. Fatty liver can promote diabetes by increasing insulin resistance.

#### Liver cirrhosis

Hepatic cirrhosis occurs as a result of permanent stress or hepatic damage. The process leading to cirrhosis can often be slowed down or even stopped. The liver has an extraordinary ability to regenerate once the underlying causes of fibrosis are controlled or treated.

#### Hepatic Carcinoma

One of the most prevalent causes of liver cancer is the infection with hepatitis B or hepatitis C virus. They are responsible for about half of all liver cancers.

## 10. Clinical Interpretation

The reference values for this assay were determined in a multicenter clinical trial. It was determined during the clinical study that:

- If the concentration of the biomarker YKL-40 (CHI3L1) in serum is <79ng/mL, this is an indication of the absence of significant fibrosis;
- if the concentration of the YKL-40 (CHI3L1) biomarker in serum is between 79-134ng/mL, this is an indication of the presence of significant fibrosis;
- if the concentration of the YKL-40 (CHI3L1) biomarker in serum is > 134ng/mL, this is an indication of the presence of liver cirrhosis;
- Possible presence of liver carcinoma and rheumatoid specimens may interfere with the determination of YKL-40 (CHI3L1) concentrations

## 11. Limitations

- The test kit is only appropriate for the detection of YKL-40 (comparable to CHI3L1) in serum samples. It is not intended for use with whole blood, heat-inactivated samples or plasma samples.
- The concentration of CHI3L1 is determined with the kit and should not be considered as absolute quantification. Hemolysis, lipemia, jaundice and microbial contamination of the samples may interfere with the correct determination of concentrations.
- Due to geographical, ethnic and age differences, each laboratory is recommended to establish its own reference value.
- The content of CHI3L1 cannot be used as the only evidence of disease.

## 12. Storage

**Storage conditions:** Store between 2 and 30°C at a dark and dry locality.

**Shelf life:** 18 months.

**Conditions of use:** After opening the aluminums packaging bag (pouch), it is strongly recommended to perform the test within 30 minutes.

**Precautions:** Do not eat, drink, or smoke while handling specimens and performing the test. All specimens are considered potentially hazardous and must be treated as infectious specimens. Follow proven precautions against microbiological risks during the test procedure. Wear protective clothing such as disposable gloves and eye protection when testing samples.

## 2. Test Accessories and Test Principle

The functioning principle: A double biomarker sandwich immunoassay is used to quantitatively determine the concentration of the biomarker YKL-40 in a serum sample. The test reagent consists mainly of the biomarker YKL-40 and the antibody labelled with colloidal gold.

The **Sample field (S)** is intended for the usage of the serum sample (100µL).

The Lateral Flow Test **Detection Strip (T)** is coated with a YKL-40 (CHI3L1) biomarker and the **Control Strip (C)** is coated with a specific polypeptide antibody (Figure 1).

After the chromatography, the immune complex is detected by the YKL-40 (CHI3L1) biomarker in the **Detection strip (T)** to form a colored red stripe and the colloidal gold-labeled polypeptide is detected by the specific polypeptide biomarker in the **Control strip (C)** to also form a colored red stripe.

Picture 1. Test Cassette with testing stripe and markings.

## 3. Sample Requirement

The assay kit is at present only usable for serum samples without microbial contamination. If a serum sample contains fibrin, particles or red blood cells, it should be centrifuged. The blood sample should be transported in accordance with national laws and regulations.

If serum samples are not used immediately, they should be stored at -20 °C. The number of freeze-thaw cycles is limited to two. When using frozen serum samples, it is important to thaw them slowly, mix well and bring to room temperature (15°-30°C). Serum samples should not be used for severe hemolysis or severe lipaemia or high bilirubin concentration.

## 13. Restrictions on Use

This **medical device** is intended for **in vitro diagnostic use under professional control**. The result is for use as a clinical reference only. Treatment of patients should be based on a combination of symptoms, clinical signs, medical history, other laboratory tests and therapeutic responses.

Before use, keep the product packaging intact. Please do not use products that exceed the validity period or are damaged.

Please follow the instructions carefully. Do not end the test halfway. If the test is half finished, it should not be resumed. This test is a disposable product. Please handle it with care to avoid danger and infection.

Due to differences in methodology or antibody specificities, assays from different manufacturers may produce different values even on exactly the same sample. These measurements are not directly comparable for clinical interpretation.

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## 4. Test Performance and Valuation

This product is a biological product. A test kit contains one test in a pouch, the sample diluent (buffer = 400 µl per bottle), these operating instructions and the calibration barcode for the reader if available (see 8. Scanner for quantitative Evaluation).

After the qualified response, the color strip on the test strip may be read by a scanner/reader to have the quantification results analyzed and displayed by the device.

## 5. Preparation of Use

1. Return the patient's test, buffer solution and specimen material to room temperature (15°-30°C).
2. Have the timer ready.
3. Do not open the sealed pouch until everything is ready to perform the test.
4. Remove the test strip seal and mark the sample ID on the test strip.
5. Place the test strip on a level table.
6. Using the transfer pipette, take 2 drops of the serum sample (100µL) and place in the application field (S).
7. Add two drops of buffer solution (100µL) and place in the application field (S). Avoid formation of air bubbles in the application field.
8. Start the stopwatch.
9. Read out the visible result exactly 15 minutes after the sample is given.

or

10. Place the cassette into the reader system and start the reader by pressing the button. The reader stops automatically by showing the quantitative result.
11. Under no circumstances should the results be displayed after 15 minutes can be read.
12. Note the results for the patient file.

### Important:

After approx. 30 seconds you can monitor the chromatography on the test strip in the cassette.

A red-pink color front forms which moves across the membrane.

In rare cases this process may be disturbed. Then please add another drop of buffer solution.

## 15. Sources/Quellen

1. Johansen, J. S., P. Christoffersen, S. Møller, P. A. Price, J. H. Henriksen, C. Garbarsch and F. Bendtsen (2000). "Serum YKL-40 is increased in patients with hepatic fibrosis." *J Hepatol* 32(6): 911-920.
2. Linlin Yan, Yongqiong Deng, Jiyuan Zhou, Hong Zhao, Guiqiang Wang (2017). "Serum YKL-40 as a biomarker for predicting significant fibrosis and advanced fibrosis in chronic hepatitis B patients with normal or mildly elevated ALT." *Hepatol Int* (2017) 11 (Suppl 1): S1-S1093.
3. Kamal, S. M., B. Turner, Q. He, J. Rasenack, L. Bianchi, A. Al Tawil, A. Nooman, M. Massoud, M. J. Koziel and N. H. Afdhal (2006). "Progression of fibrosis in hepatitis C with and without schistosomiasis: correlation with serum markers of fibrosis." *Hepatology* 43(4): 771-779.
4. Wang J, Gao F, Mo F, Hong X, Wang H, Zheng S, Lin B (2009). "Identification of for Liver cancer through integrative secretome and transcriptome analysis." *Proteomics-Clinical Applications*, 3: 541-551, 2009.
5. Hood, L., J. R. Heath, M. E. Phelps and B. Lin (2004). "Systems biology and new technologies enable predictive and preventative medicine." *Science* 306(5696): 640-643.
6. Huang, H., T. Wu, J. Mao, Y. Fang, J. Zhang, L. Wu, S. Zheng, B. Lin and H. Pan (2015). "CHI3L1 Is a Liver-Enriched, Noninvasive Biomarker That Can Be Used to Stage and Diagnose Substantial Hepatic Fibrosis." *OMICS, a Journal of Integrative Biology* 19(6): 339-345.
7. Kumagai, E., Y. Mano, S. Yoshio, H. Shoji, M. Sugiyama, M. Korenaga, T. Ishida, T. Arai, N. Itokawa, M. Atsukawa, H. Hyogo, K. Chayama, T. Ohashi, K. Ito, M. Yoneda, T. Kawaguchi, T. Torimura, Y. Nozaki, S. Watanabe, M. Mizokami and T. Kanto (2016). "Serum YKL-40 as a marker of liver fibrosis in patients with non-alcoholic fatty liver disease." *Sci Rep* 6: 35282.
8. Lin Wang, Yameng Sun, Hong You, Jidong Jia (2017). "The decrease of serum CHI3L1 levels after antiviral therapy can monitor reversal of liver fibrosis." *Hepatol Int* (2017) 11 (Suppl 1): S1-S1093.
9. Hakala, B. E., C. White and A. D. Recklies (1993). "Human cartilage gp-39, a major secretory product of articular chondrocytes and synovial cells, is a mammalian member of a chitinase protein family." *J Biol Chem* 268(34): 25803-25810.

Please visit us at [www.PharmACT.eu](http://www.PharmACT.eu)